



TELO2 rabbit pAb

Cat No.:ES19571

For research use only

Overview

Product Name	TELO2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TELO2 AA range: 605-655
Specificity	This antibody detects endogenous levels of human TELO2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Telomere length regulation protein TEL2 homolog (Protein clk-2 homolog) (hCLK2)
Gene Name	TELO2 KIAA0683
Cellular localization	Cytoplasm. Membrane. Nucleus. Chromosome, telomere .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9894
Human Swiss-Prot Number	Q9Y4R8
Alternative Names	
Background	

