

# DOCK9 rabbit pAb

Cat No.:ES19573

For research use only

## Overview

Product Name	DOCK9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DOCK9 AA range: 1516-1566
Specificity	This antibody detects endogenous levels of human DOCK9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Dedicator of cytokinesis protein 9 (Cdc42 guanine nucleotide exchange factor zizimin-1)
Gene Name	DOCK9 KIAA1058
Cellular localization	Endomembrane system . Associated with membranes. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23348
Human Swiss-Prot Number	Q9BZ29
Alternative Names	
Background	

