

MACD2 rabbit pAb

Cat No.:ES19587

For research use only

Overview

Product Name	MACD2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MACD2 AA range: 166-216
Specificity	This antibody detects endogenous levels of human MACD2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	O-acetyl-ADP-ribose deacetylase MACROD2 (EC 3.5.1.-) (MACRO domain-containing protein 2)
Gene Name	MACROD2 C20orf133
Cellular localization	Nucleus . Recruited to DNA lesions, probably via mono-APD-ribosylated proteins. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	140733
Human Swiss-Prot Number	A1Z1Q3
Alternative Names	
Background	

