



HEAT4 rabbit pAb

Cat No.:ES19588

For research use only

Overview

Product Name	HEAT4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HEAT4 AA range: 502-552
Specificity	This antibody detects endogenous levels of human HEAT4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	HEAT repeat-containing protein 4
Gene Name	HEATR4
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	399671
Human Swiss-Prot Number	Q86WZ0
Alternative Names	
Background	

