



QPCTL rabbit pAb

Cat No.:ES19590

For research use only

Overview

Product Name	QPCTL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human QPCTL AA range: 152-202
Specificity	This antibody detects endogenous levels of human QPCTL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Glutaminy-peptide cyclotransferase-like protein (EC 2.3.2.5) (Golgi-resident glutaminy-peptide cyclotransferase) (isoQC) (gQC)
Gene Name	QPCTL
Cellular localization	Golgi apparatus membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54814
Human Swiss-Prot Number	Q9NXS2
Alternative Names	
Background	

