



MOS1 rabbit pAb

Cat No.:ES19594

For research use only

Overview

Product Name	MOS1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MOS1 AA range: 9-59
Specificity	This antibody detects endogenous levels of human MOS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Mitochondrial inner membrane organizing system protein 1
Gene Name	MINOS1 C1orf151
Cellular localization	Mitochondrion inner membrane ; Single-pass membrane protein . The C-terminus is located in the intermembrane space (By similarity), while the location of the N-terminus has not been determined yet. As some programs predict the presence of 2 closely apposed membrane domains, it has been proposed that the protein may cross the membrane twice and that both termini may face the intermembrane space (PubMed:22114354). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	440574
Human Swiss-Prot Number	Q5TGZ0
Alternative Names	





ELK Biotechnology

Background



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C