

NFIA rabbit pAb

Cat No.:ES19598

For research use only

Overview

Product Name	NFIA rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NFIA AA range: 391-441
Specificity	This antibody detects endogenous levels of human NFIA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nuclear factor 1 A-type (NF1-A) (Nuclear factor 1/A) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/A) (NF-I/A) (NFI-A) (TGGCA-binding protein)
Gene Name	NFIA KIAA1439
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4774
Human Swiss-Prot Number	Q12857
Alternative Names	
Background	

