

SHSA9 rabbit pAb

Cat No.:ES19601

For research use only

Overview

Product Name	SHSA9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SHSA9 AA range: 352-402
Specificity	This antibody detects endogenous levels of human SHSA9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein shisa-9
Gene Name	SHISA9
Cellular localization	Cell projection, dendritic spine membrane ; Single-pass type I membrane protein . Cell junction, synapse .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	729993
Human Swiss-Prot Number	B4DS77
Alternative Names	
Background	

