



HXD4 rabbit pAb

Cat No.:ES19602

For research use only

Overview

Product Name	HXD4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HXD4 AA range: 99-149
Specificity	This antibody detects endogenous levels of human HXD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein Hox-D4 (Homeobox protein HHO.C13) (Homeobox protein Hox-4B) (Homeobox protein Hox-5.1)
Gene Name	HOXD4 HOX4B
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3233
Human Swiss-Prot Number	P09016
Alternative Names	
Background	

