

MF2L2 rabbit pAb

Cat No.:ES19613

For research use only

Overview

Product Name	MF2L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MF2L2 AA range: 97-147
Specificity	This antibody detects endogenous levels of human MF2L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Probable guanine nucleotide exchange factor MCF2L2 (Dbs-related Rho family guanine nucleotide exchange factor) (MCF2-transforming sequence-like protein 2)
Gene Name	MCF2L2 DRG KIAA0861
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23101
Human Swiss-Prot Number	Q86YR7
Alternative Names	
Background	

