

PERQ1 rabbit pAb

Cat No.:ES19616

For research use only

Overview

Product Name	PERQ1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PERQ1 AA range: 684-734
Specificity	This antibody detects endogenous levels of human PERQ1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	PERQ amino acid-rich with GYF domain-containing protein 1 (GRB10-interacting GYF protein 1)
Gene Name	GIGYF1 CDS2 PERQ1 PP3360
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	64599
Human Swiss-Prot Number	O75420
Alternative Names	
Background	

