



ALKB6 rabbit pAb

Cat No.:ES19617

For research use only

Overview

Product Name	ALKB6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ALKB6 AA range: 17-67
Specificity	This antibody detects endogenous levels of human ALKB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Probable alpha-ketoglutarate-dependent dioxygenase ABH6 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 6)
Gene Name	ALKBH6 ABH6
Cellular localization	Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84964
Human Swiss-Prot Number	Q3KRA9
Alternative Names	
Background	

