

RN19B rabbit pAb

Cat No.:ES19622

For research use only

Overview

Product Name	RN19B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RN19B AA range: 272-322
Specificity	This antibody detects endogenous levels of human RN19B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase RNF19B (EC 6.3.2.-) (IBR domain-containing protein 3) (Natural killer lytic-associated molecule) (RING finger protein 19B)
Gene Name	RNF19B IBRD3 NKLAM
Cellular localization	Cytoplasmic granule membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	127544
Human Swiss-Prot Number	Q6ZMZ0
Alternative Names	
Background	

