



# CERS6 rabbit pAb

Cat No.:ES19625

For research use only

## Overview

Product Name	CERS6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CERS6 AA range: 73-123
Specificity	This antibody detects endogenous levels of human CERS6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ceramide synthase 6 (CerS5) (LAG1 longevity assurance homolog 6)
Gene Name	CERS6 LASS6
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	253782
Human Swiss-Prot Number	Q6ZMG9
Alternative Names	
Background	

