

# NBPFF rabbit pAb

Cat No.:ES19626

For research use only

## Overview

Product Name	NBPFF rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NBPFF AA range: 382-432
Specificity	This antibody detects endogenous levels of human NBPFF
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Neuroblastoma breakpoint family member 15
Gene Name	NBPF15
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	284565
Human Swiss-Prot Number	Q8N660
Alternative Names	
Background	

