

# LCE3B rabbit pAb

Cat No.:ES19627

For research use only

## Overview

Product Name	LCE3B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LCE3B AA range: 8-58
Specificity	This antibody detects endogenous levels of human LCE3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Late cornified envelope protein 3B (Late envelope protein 14)
Gene Name	LCE3B LEP14
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	353143
Human Swiss-Prot Number	Q5TA77
Alternative Names	
Background	

