

TMHS rabbit pAb

Cat No.:ES19631

For research use only

Overview

Product Name	TMHS rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TMHS AA range: 28-78
Specificity	This antibody detects endogenous levels of human TMHS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tetraspan membrane protein of hair cell stereocilia (Lipoma HMGIC fusion partner-like 5 protein)
Gene Name	LHFPL5 TMHS
Cellular localization	Cell membrane ; Multi-pass membrane protein . Efficient localization to the plasma membrane requires the presence of PCDH15. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	222662
Human Swiss-Prot Number	Q8TAF8
Alternative Names	
Background	

