



LMBL3 rabbit pAb

Cat No.:ES19634

For research use only

Overview

Product Name	LMBL3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LMBL3 AA range: 128-178
Specificity	This antibody detects endogenous levels of human LMBL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Lethal(3)malignant brain tumor-like protein 3 (H-l(3)mbt-like protein 3) (L(3)mbt-like protein 3) (MBT-1)
Gene Name	L3MBTL3 KIAA1798 MBT1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84456
Human Swiss-Prot Number	Q96JM7
Alternative Names	
Background	

