



YRDC rabbit pAb

Cat No.:ES19641

For research use only

Overview

Product Name	YRDC rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human YRDC AA range: 214-264
Specificity	This antibody detects endogenous levels of human YRDC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	YrdC domain-containing protein, mitochondrial (Dopamine receptor-interacting protein 3) (Ischemia/reperfusion-inducible protein homolog) (hIRIP)
Gene Name	YRDC DRIP3 IRIP
Cellular localization	Cytoplasm . Mitochondrion . Cell membrane ; Peripheral membrane protein . A large fraction localizes in the cytoplasm, whereas a smaller fraction is imported to mitochondria. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79693
Human Swiss-Prot Number	Q86U90
Alternative Names	
Background	

