

FBW10 rabbit pAb

Cat No.:ES19643

For research use only

Overview

Product Name	FBW10 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FBW10 AA range: 779-829
Specificity	This antibody detects endogenous levels of human FBW10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	F-box/WD repeat-containing protein 10 (F-box and WD-40 domain-containing protein 10) (Ubiquitin ligase-specificity factor)
Gene Name	FBXW10
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10517
Human Swiss-Prot Number	Q5XX13
Alternative Names	
Background	

