



FCF1 rabbit pAb

Cat No.:ES19644

For research use only

Overview

Product Name	FCF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FCF1 AA range: 56-106
Specificity	This antibody detects endogenous levels of human FCF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	rRNA-processing protein FCF1 homolog
Gene Name	FCF1 C14orf111 CGI-35
Cellular localization	Nucleus, nucleolus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51077
Human Swiss-Prot Number	Q9Y324
Alternative Names	
Background	

