



NAF1 rabbit pAb

Cat No.:ES19646

For research use only

Overview

Product Name	NAF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NAF1 AA range: 409-459
Specificity	This antibody detects endogenous levels of human NAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	H/ACA ribonucleoprotein complex non-core subunit NAF1 (hNAF1)
Gene Name	NAF1
Cellular localization	Cytoplasm. Nucleus. Shuttles between the cytoplasm and the nucleus. Absent from the nucleolus (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	92345
Human Swiss-Prot Number	Q96HR8
Alternative Names	
Background	

