



TGM7 rabbit pAb

Cat No.:ES19647

For research use only

Overview

Product Name	TGM7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TGM7 AA range: 112-162
Specificity	This antibody detects endogenous levels of human TGM7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Protein-glutamine gamma-glutamyltransferase Z (Transglutaminase Z) (TG(Z)) (TGZ) (TGase Z) (EC 2.3.2.13) (Transglutaminase-7) (TGase-7)
Gene Name	TGM7
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	116179
Human Swiss-Prot Number	Q96PF1
Alternative Names	
Background	

