



# C1AS1 rabbit pAb

Cat No.:ES19654

For research use only

## Overview

Product Name	C1AS1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human C1AS1 AA range: 19-69
Specificity	This antibody detects endogenous levels of human C1AS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative uncharacterized protein C17orf76-AS1, mitochondrial (C17orf76 antisense RNA 1) (C17orf76 antisense gene protein 1)
Gene Name	C17orf76-AS1 C17orf45 NCRNA00188
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q8N1F1
Alternative Names	
Background	

