



S15A2 rabbit pAb

Cat No.:ES19655

For research use only

Overview

Product Name	S15A2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human S15A2 AA range: 484-534
Specificity	This antibody detects endogenous levels of human S15A2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Solute carrier family 15 member 2 (Kidney H(+)/peptide cotransporter) (Oligopeptide transporter, kidney isoform) (Peptide transporter 2)
Gene Name	SLC15A2 PEPT2
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6565
Human Swiss-Prot Number	Q16348
Alternative Names	
Background	

