

# BABA1 rabbit pAb

Cat No.:ES19657

For research use only

## Overview

Product Name	BABA1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human BABA1 AA range: 84-134
Specificity	This antibody detects endogenous levels of human BABA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	BRISC and BRCA1-A complex member 1 (Mediator of RAP80 interactions and targeting subunit of 40 kDa) (New component of the BRCA1-A complex)
Gene Name	BABAM1 C19orf62 MERIT40 NBA1 HSPC142
Cellular localization	Cytoplasm . Nucleus . Localizes at sites of DNA damage at double-strand breaks (DSBs). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29086
Human Swiss-Prot Number	Q9NWX8
Alternative Names	
Background	

