

ZMYM1 rabbit pAb

Cat No.:ES19658

For research use only

Overview

Product Name	ZMYM1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZMYM1 AA range: 904-954
Specificity	This antibody detects endogenous levels of human ZMYM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger MYM-type protein 1
Gene Name	ZMYM1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79830
Human Swiss-Prot Number	Q5SVZ6
Alternative Names	
Background	

