

DB127 rabbit pAb

Cat No.:ES19659

For research use only

Overview

Product Name	DB127 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DB127 AA range: 44-94
Specificity	This antibody detects endogenous levels of human DB127
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Beta-defensin 127 (Beta-defensin 27) (DEFB-27) (Defensin, beta 127)
Gene Name	DEFB127 C20orf73 DEFB27 UNQ1956/PRO6071
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	140850
Human Swiss-Prot Number	Q9H1M4
Alternative Names	
Background	

