

CC062 rabbit pAb

Cat No.:ES19660

For research use only

Overview

Product Name	CC062 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CC062 AA range: 173-223
Specificity	This antibody detects endogenous levels of human CC062
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C3orf62
Gene Name	C3orf62
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	375341
Human Swiss-Prot Number	Q6ZUJ4
Alternative Names	
Background	

