



# PRCD rabbit pAb

Cat No.:ES19661

For research use only

## Overview

Product Name	PRCD rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PRCD AA range: 4-54
Specificity	This antibody detects endogenous levels of human PRCD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Progressive rod-cone degeneration protein
Gene Name	PRCD
Cellular localization	Cell projection, cilium, photoreceptor outer segment . Membrane ; Lipid-anchor ; Cytoplasmic side . Endoplasmic reticulum . Golgi apparatus . Localizes to photoreceptor disk membranes in the photoreceptor outer segment (PubMed:27613864). The secretion in media described in PubMed:24992209 is probably an experimental artifact (PubMed:24992209). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	768206
Human Swiss-Prot Number	Q00LT1
Alternative Names	
Background	

