

DGCR6 rabbit pAb

Cat No.:ES19667

For research use only

Overview

Product Name	DGCR6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DGCR6 AA range: 34-84
Specificity	This antibody detects endogenous levels of human DGCR6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein DGCR6 (DiGeorge syndrome critical region 6)
Gene Name	DGCR6
Cellular localization	Nucleus . Predominantly nuclear.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8214
Human Swiss-Prot Number	Q14129
Alternative Names	
Background	

