



VAT1L rabbit pAb

Cat No.:ES19674

For research use only

Overview

Product Name	VAT1L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VAT1L AA range: 118-168
Specificity	This antibody detects endogenous levels of human VAT1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Synaptic vesicle membrane protein VAT-1 homolog-like (EC 1.-.-.-)
Gene Name	VAT1L KIAA1576
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57687
Human Swiss-Prot Number	Q9HCJ6
Alternative Names	
Background	

