



# MS4A5 rabbit pAb

Cat No.:ES19675

For research use only

## Overview

<b>Product Name</b>	MS4A5 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human MS4A5 AA range: 46-96
<b>Specificity</b>	This antibody detects endogenous levels of human MS4A5
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Membrane-spanning 4-domains subfamily A member 5 (CD20 antigen-like 2) (Testis-expressed transmembrane protein 4)
<b>Gene Name</b>	MS4A5 CD20L2 TETM4
<b>Cellular localization</b>	Membrane; Multi-pass membrane protein.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	64232
<b>Human Swiss-Prot Number</b>	Q9H3V2
<b>Alternative Names</b>	
<b>Background</b>	

