

# F131C rabbit pAb

Cat No.:ES19676

For research use only

## Overview

Product Name	F131C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human F131C AA range: 130-180
Specificity	This antibody detects endogenous levels of human F131C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM131C
Gene Name	FAM131C C1orf117
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	348487
Human Swiss-Prot Number	Q96AQ9
Alternative Names	
Background	

