

DUS8 rabbit pAb

Cat No.:ES19678

For research use only

Overview

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name | DUS8 rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human DUS8 AA range: 98-148 |
| Specificity | This antibody detects endogenous levels of human DUS8 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Dual specificity protein phosphatase 8 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH-5) |
| Gene Name | DUSP8 C11orf81 VH5 |
| Cellular localization | Cytoplasm . Nucleus . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 1850 |
| Human Swiss-Prot Number | Q13202 |
| Alternative Names | |
| Background | |

