



ZSC21 rabbit pAb

Cat No.:ES19679

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | ZSC21 rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human ZSC21 AA range: 222-272 |
| Specificity | This antibody detects endogenous levels of human ZSC21 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Zinc finger and SCAN domain-containing protein 21 (Renal carcinoma antigen NY-REN-21) (Zinc finger protein 38 homolog) (Zfp-38) |
| Gene Name | ZSCAN21 ZFP38 ZNF38 |
| Cellular localization | Nucleus . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 7589 |
| Human Swiss-Prot Number | Q9Y5A6 |
| Alternative Names | |
| Background | |

