



TPK1 rabbit pAb

Cat No.:ES19681

For research use only

Overview

Product Name	TPK1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TPK1 AA range: 171-221
Specificity	This antibody detects endogenous levels of human TPK1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Thiamin pyrophosphokinase 1 (hTPK1) (EC 2.7.6.2) (Placental protein 20) (PP20) (Thiamine pyrophosphokinase 1)
Gene Name	TPK1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	27010
Human Swiss-Prot Number	Q9H3S4
Alternative Names	
Background	

