



P4R3B rabbit pAb

Cat No.:ES19687

For research use only

Overview

Product Name	P4R3B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human P4R3B AA range: 425-475
Specificity	This antibody detects endogenous levels of human P4R3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Serine/threonine-protein phosphatase 4 regulatory subunit 3B (SMEK homolog 2)
Gene Name	SMEK2 KIAA1387 PP4R3B PPP4R3B
Cellular localization	Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Nucleus . In interphase localized in the cytoplasm and (with higher levels) the nucleus. During metaphase located in pericentriolar regions.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57223
Human Swiss-Prot Number	Q5MIZ7
Alternative Names	
Background	

