



VOPP1 rabbit pAb

Cat No.:ES19696

For research use only

Overview

Product Name	VOPP1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human VOPP1 AA range: 120-170
Specificity	This antibody detects endogenous levels of human VOPP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Vesicular, overexpressed in cancer, prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N)
Gene Name	VOPP1 ECOP GASP
Cellular localization	Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . When overexpressed, localizes in the nucleus and perinuclear regions.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	81552
Human Swiss-Prot Number	Q96AW1
Alternative Names	
Background	

