



ARL3 rabbit pAb

Cat No.:ES19698

For research use only

Overview

Product Name	ARL3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ARL3 AA range: 109-159
Specificity	This antibody detects endogenous levels of human ARL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	ADP-ribosylation factor-like protein 3
Gene Name	ARL3 ARFL3
Cellular localization	Golgi apparatus membrane; Peripheral membrane protein; Cytoplasmic side. Cytoplasm, cytoskeleton, spindle. Nucleus. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm. Cell projection, cilium . Detected predominantly in the photoreceptor connecting cilium. Present on the mitotic spindle. Centrosome-associated throughout the cell cycle. Not detected to interphase microtubules.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	403
Human Swiss-Prot Number	P36405
Alternative Names	
Background	

