

LBR rabbit pAb

Cat No.:ES19700

For research use only

Overview

Product Name	LBR rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LBR AA range: 411-461
Specificity	This antibody detects endogenous levels of human LBR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lamin-B receptor (Integral nuclear envelope inner membrane protein) (LMN2R)
Gene Name	LBR
Cellular localization	Nucleus inner membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane . Cytoplasm . Nucleus . Nucleus; nuclear rim. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3930
Human Swiss-Prot Number	Q14739
Alternative Names	
Background	

