

ST1E1 rabbit pAb

Cat No.:ES19704

For research use only

Overview

Product Name	ST1E1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ST1E1 AA range: 209-259
Specificity	This antibody detects endogenous levels of human ST1E1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Estrogen sulfotransferase (EC 2.8.2.4) (EST-1) (Sulfotransferase 1E1) (ST1E1) (Sulfotransferase, estrogen-preferring)
Gene Name	SULT1E1 STE
Cellular localization	Cytoplasm, cytosol .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6783
Human Swiss-Prot Number	P49888
Alternative Names	
Background	

