



ZN746 rabbit pAb

Cat No.:ES19708

For research use only

Overview

Product Name	ZN746 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZN746 AA range: 37-87
Specificity	This antibody detects endogenous levels of human ZN746
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 746 (Parkin-interacting substrate) (PARIS)
Gene Name	ZNF746 PARIS
Cellular localization	Cytoplasm . Nucleus . Mainly localizes to the cytoplasm; probably translocates to the nucleus to repress selected genes.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	155061
Human Swiss-Prot Number	Q6NUN9
Alternative Names	
Background	

