

PGTB2 rabbit pAb

Cat No.:ES19710

For research use only

Overview

Product Name	PGTB2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PGTB2 AA range: 64-114
Specificity	This antibody detects endogenous levels of human PGTB2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Geranylgeranyl transferase type-2 subunit beta (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit beta) (GGTase-II-beta) (Rab geranyl-geranyltransferase subunit beta) (Rab GG transferase beta)
Gene Name	RABGGTB GGTB
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5876
Human Swiss-Prot Number	P53611
Alternative Names	
Background	

