



TPC10 rabbit pAb

Cat No.:ES19711

For research use only

Overview

Product Name	TPC10 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TPC10 AA range: 468-518
Specificity	This antibody detects endogenous levels of human TPC10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Trafficking protein particle complex subunit 10 (Epilepsy holoprosencephaly candidate 1 protein) (EHOC-1) (Protein GT334) (Trafficking protein particle complex subunit TMEM1) (Transport protein partic
Gene Name	TRAPPC10 EHOC1 TMEM1
Cellular localization	Golgi apparatus, cis-Golgi network .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7109
Human Swiss-Prot Number	P48553
Alternative Names	
Background	

