



ZN214 rabbit pAb

Cat No.:ES19712

For research use only

Overview

Product Name	ZN214 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZN214 AA range: 509-559
Specificity	This antibody detects endogenous levels of human ZN214
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 214 (BWSCR2-associated zinc finger protein 1) (BAZ-1)
Gene Name	ZNF214 BAZ1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7761
Human Swiss-Prot Number	Q9UL59
Alternative Names	
Background	

