



GLOD4 rabbit pAb

Cat No.:ES19714

For research use only

Overview

Product Name	GLOD4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GLOD4 AA range: 53-103
Specificity	This antibody detects endogenous levels of human GLOD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Glyoxalase domain-containing protein 4
Gene Name	GLOD4 C17orf25 CGI-150 My027
Cellular localization	Mitochondrion .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51031
Human Swiss-Prot Number	Q9HC38
Alternative Names	
Background	

