



ZF64B rabbit pAb

Cat No.:ES19717

For research use only

Overview

Product Name	ZF64B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZF64B AA range: 292-342
Specificity	This antibody detects endogenous levels of human ZF64B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 64 homolog, isoforms 3 and 4 (Zfp-64) (Zinc finger protein 338)
Gene Name	ZFP64 ZNF338
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55734
Human Swiss-Prot Number	Q9NTW7
Alternative Names	
Background	

