

# FA58B rabbit pAb

Cat No.:ES19721

For research use only

## Overview

Product Name	FA58B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FA58B AA range: 74-124
Specificity	This antibody detects endogenous levels of human FA58B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	Putative cyclin-related protein FAM58B
Gene Name	FAM58BP FAM58B
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	339521
Human Swiss-Prot Number	P0C7Q3
Alternative Names	
Background	

