

# LY66C rabbit pAb

Cat No.:ES19723

For research use only

## Overview

Product Name	LY66C rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LY66C AA range: 53-103
Specificity	This antibody detects endogenous levels of human LY66C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lymphocyte antigen 6 complex locus protein G6c
Gene Name	LY6G6C C6orf24 G6C NG24 UNQ1947/PRO4430
Cellular localization	Cell membrane ; Lipid-anchor, GPI-anchor .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	80740
Human Swiss-Prot Number	O95867
Alternative Names	
Background	

